

## Bag-3 Polyclonal Antibody

Catalog No: #41910

Package Size: #41910-1 50ul #41910-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bag-3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	Bag-3 Polyclonal Antibody detects endogenous levels of Bag-3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Bag-3.
Target Name	Bag-3
Other Names	BAG3; BIS; BAG family molecular chaperone regulator 3; BAG-3; Bcl-2-associated athanogene 3; Bcl-2-binding protein Bis; Docking protein CAIR-1
Accession No.	Swiss-Prot: O95817NCBI Gene ID: 9531
Uniprot	O95817
GeneID	9531;
SDS-PAGE MW	65kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

## Images



Western Blot analysis of A549 cells using Bag-3 Polyclonal Antibody.

---

Note: This product is for in vitro research use only