

Tenascin-C Polyclonal Antibody

Catalog No: #41914

Package Size: #41914-1 50ul #41914-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tenascin-C Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Tenascin-C Polyclonal Antibody detects endogenous levels of Tenascin-C protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Tenascin-C.
Target Name	Tenascin-C
Other Names	TNC; HXB; Tenascin; TN; Cytotactin; GMEM; GP 150-225; Glioma-associated-extracellular matrix antigen; Hexabrachion; JI; Myotendinous antigen; Neuronectin; Tenascin-C; TN-C
Accession No.	Swiss-Prot: P24821NCBI Gene ID: 3371
Uniprot	P24821
GeneID	3371;
SDS-PAGE MW	240kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of KB cells using Tenascin-C Polyclonal Antibody

Note: This product is for in vitro research use only