

PSA Polyclonal Antibody

Catalog No: #41917

Package Size: #41917-1 50ul #41917-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PSA Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF
Species Reactivity	Human Rat
Specificity	PSA Polyclonal Antibody detects endogenous levels of PSA protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human PSA.
Target Name	PSA
Other Names	KLK3; APS; Prostate-specific antigen; PSA; Gamma-seminoprotein; Semin; Kallikrein-3; P-30 antigen; Semenogelase
Accession No.	Swiss-Prot: P07288NCBI Gene ID: 354
Uniprot	P07288
GeneID	354;
SDS-PAGE MW	34kDa
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500-1:2000


IHC 1:100-1:300

IF 1:200-1:1000

Images



Western blot analysis of lysates from mouse heart cells, using Prostate-specific Antigen Antibody. The lane on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human prostate carcinoma tissue, using Prostate-specific Antigen Antibody. The picture on the right is blocked with the synthesized peptide.

Note: This product is for in vitro research use only