

FGF-17 Polyclonal Antibody

Catalog No: #41920

Package Size: #41920-1 50ul #41920-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FGF-17 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	FGF-17 Polyclonal Antibody detects endogenous levels of FGF-17 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human FGF-17.
Target Name	FGF-17
Other Names	FGF17; Fibroblast growth factor 17; FGF-17
Accession No.	Swiss-Prot: O60258NCBI Gene ID: 8822
Uniprot	O60258
GeneID	8822;
SDS-PAGE MW	25kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of NIH-3T3 cells using FGF-17 Polyclonal Antibody

Note: This product is for in vitro research use only