

NOXIN Polyclonal Antibody

Catalog No: #41930

Package Size: #41930-1 50ul #41930-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NOXIN Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	NOXIN Polyclonal Antibody detects endogenous levels of NOXIN protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human NOXIN.
Target Name	NOXIN
Other Names	NOXIN; C11orf82; Nitric oxide-inducible gene protein
Accession No.	Swiss-Prot: Q8IXT1NCBI Gene ID: 220042
Uniprot	Q8IXT1
GeneID	220042;
SDS-PAGE MW	85kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of 823 293T AD293 Hela HepG2 cells using NOXIN Polyclonal Antibody

Note: This product is for in vitro research use only