

IGFBP1 Polyclonal Antibody

Catalog No: #41933

Package Size: #41933-1 50ul #41933-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IGFBP1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	IGFBP1 Polyclonal Antibody detects endogenous levels of IGFBP1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human IGFBP1.
Target Name	IGFBP1
Other Names	IGFBP1; IBP1; Insulin-like growth factor-binding protein 1; IBP-1; IGF-binding protein 1; IGFBP-1; Placental protein 12; PP12
Accession No.	Swiss-Prot: P08833NCBI Gene ID: 3484
Uniprot	P08833
GeneID	3484;
SDS-PAGE MW	28kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of MCF7 cells using IGFBP1 Polyclonal Antibody.

Note: This product is for in vitro research use only