

CD89 Polyclonal Antibody

Catalog No: #41934

Package Size: #41934-1 50ul #41934-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD89 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	CD89 Polyclonal Antibody detects endogenous levels of CD89 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CD89.
Target Name	CD89
Other Names	FCAR; CD89; Immunoglobulin alpha Fc receptor; IgA Fc receptor; CD89
Accession No.	Swiss-Prot: P24071NCBI Gene ID: 2204
Uniprot	P24071
GeneID	2204;
SDS-PAGE MW	32kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of NIH-3T3 cells using CD89 Polyclonal Antibody

Note: This product is for in vitro research use only