

GABAA R α 1 Polyclonal Antibody

Catalog No: #41935

Package Size: #41935-1 50ul #41935-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GABAA R α 1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA IF
Species Reactivity	Hu Ms Rt
Specificity	GABAA R α 1 Polyclonal Antibody detects endogenous levels of GABAA R α 1 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human GABAA R α 1.
Target Name	GABAA R α 1
Other Names	GABRA1; Gamma-aminobutyric acid receptor subunit alpha-1; GABA(A) receptor subunit alpha-1
Accession No.	Swiss-Prot: P14867NCBI Gene ID: 2554
Uniprot	P14867
GeneID	2554;
SDS-PAGE MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/10000.


Not yet tested in other applications.

IF: 1:50-200

IHC: 1:100-1:300

Images

Western Blot analysis of NIH-3T3 KB cells using GABAA R α 1 Polyclonal Antibody



Immunohistochemical analysis of paraffin-embedded rat-brain, antibody was diluted at 1:100

Note: This product is for in vitro research use only