

## PSGL-1 Polyclonal Antibody

Catalog No: #41938

Package Size: #41938-1 50ul #41938-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	PSGL-1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	PSGL-1 Polyclonal Antibody detects endogenous levels of PSGL-1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human PSGL-1.
Target Name	PSGL-1
Other Names	SELPLG; P-selectin glycoprotein ligand 1; PSGL-1; Selectin P ligand; CD162
Accession No.	Swiss-Prot: Q14242NCBI Gene ID: 6404
Uniprot	Q14242
GeneID	6404;
SDS-PAGE MW	45kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

## Images



Western Blot analysis of NIH-3T3 cells using PSGL-1 Polyclonal Antibody

---

Note: This product is for in vitro research use only