

CD53 Polyclonal Antibody

Catalog No: #41956

Package Size: #41956-1 50ul #41956-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD53 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	CD53 Polyclonal Antibody detects endogenous levels of CD53 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CD53.
Target Name	CD53
Other Names	CD53; MOX44; TSPAN25; Leukocyte surface antigen CD53; Cell surface glycoprotein CD53; Tetraspanin-25; Tspan-25; CD53
Accession No.	Swiss-Prot: P19397NCBI Gene ID: 963
Uniprot	P19397
GeneID	963;
SDS-PAGE MW	24kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of KB cells using CD53 Polyclonal Antibody



Western Blot analysis of KB cells using CD53 Polyclonal Antibody

Note: This product is for in vitro research use only