

# Peroxin 13 antibody

Catalog No: #22718

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	Peroxin 13 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 342 and 403 of Human PEX13
Target Name	Peroxin 13
Accession No.	Swiss-Prot:Q92968Gene ID:5194
Uniprot	Q92968
GeneID	5194;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details

Predicted MW: 44kd

Western blotting: 1:500-1:3000

## Images



Sample(30 ug of whole cell lysate)  
A: MOLT4  
7.5% SDS PAGE  
Primary antibody diluted at 1: 1000

## Background

This gene encodes a peroxisomal membrane protein that binds the type 1 peroxisomal targeting signal receptor via a SH3 domain located in the cytoplasm. Mutations and deficiencies in peroxisomal protein importing and peroxisome assembly lead to peroxisomal biogenesis disorders, an

example of which is Zellweger syndrome. [provided by RefSeq]

---

Note: This product is for in vitro research use only