

Integrin $\beta 3$ Polyclonal Antibody

Catalog No: #41957

Package Size: #41957-1 50ul #41957-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Integrin $\beta 3$ Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	Integrin $\beta 3$ Polyclonal Antibody detects endogenous levels of Integrin $\beta 3$ protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Integrin $\beta 3$.
Target Name	Integrin $\beta 3$
Other Names	ITGB3; GP3A; Integrin beta-3; Platelet membrane glycoprotein IIIa; GPIIIa; CD61
Accession No.	Swiss-Prot: P05106NCBI Gene ID: 3690
Uniprot	P05106
GeneID	3690;
SDS-PAGE MW	87kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details


Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images

Western Blot analysis of 3T3 cells using Integrin $\beta 3$ Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using Integrin β 3 Polyclonal Antibody

Note: This product is for in vitro research use only