

CD85d Polyclonal Antibody

Catalog No: #41960



Package Size: #41960-1 50ul #41960-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD85d Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	CD85d Polyclonal Antibody detects endogenous levels of CD85d protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CD85d.
Target Name	CD85d
Other Names	LILRB2; ILT4; LIR2; MIR10; Leukocyte immunoglobulin-like receptor subfamily B member 2; LIR-2; Leukocyte immunoglobulin-like receptor 2; CD85 antigen-like family member D; Immunoglobulin-like transcript 4; ILT-4; Monocyte/macrophage immunoglobulin-like re
Accession No.	Swiss-Prot: Q8N423NCBI Gene ID: 10288
Uniprot	Q8N423
GeneID	10288;
SDS-PAGE MW	66kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using CD85d Polyclonal Antibody

Note: This product is for in vitro research use only