

Nectin 2 Polyclonal Antibody

Catalog No: #41968



Package Size: #41968-1 50ul #41968-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Nectin 2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	Nectin 2 Polyclonal Antibody detects endogenous levels of Nectin 2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Nectin 2.
Target Name	Nectin 2
Other Names	PVRL2; HVEB; PRR2; Poliovirus receptor-related protein 2; Herpes virus entry mediator B; Herpesvirus entry mediator B; HveB; Nectin-2; CD112
Accession No.	Swiss-Prot: Q92692NCBI Gene ID: 5819
Uniprot	Q92692
GeneID	5819;
SDS-PAGE MW	65kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.


Not yet tested in other applications.

Images



Western Blot analysis of MOUSE-BRAIN cells using Nectin 2 Polyclonal Antibody

Western Blot analysis of mouse brain cells using Nectin 2 Polyclonal Antibody



Note: This product is for in vitro research use only