

## LIFR Polyclonal Antibody

Catalog No: #41969

Package Size: #41969-1 50ul #41969-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	LIFR Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	LIFR Polyclonal Antibody detects endogenous levels of LIFR protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human LIFR.
Target Name	LIFR
Other Names	LIFR; Leukemia inhibitory factor receptor; LIF receptor; LIF-R; CD118
Accession No.	Swiss-Prot: P42702NCBI Gene ID: 3977
Uniprot	P42702
GeneID	3977;
SDS-PAGE MW	120kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

## Images



Western Blot analysis of 293 cells using LIFR Polyclonal Antibody

---

Note: This product is for in vitro research use only