

BST-2 Polyclonal Antibody

Catalog No: #41977

Package Size: #41977-1 50ul #41977-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	BST-2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	BST-2 Polyclonal Antibody detects endogenous levels of BST-2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human BST-2.
Target Name	BST-2
Other Names	BST2; Bone marrow stromal antigen 2; BST-2; HM1.24 antigen; Tetherin; CD317
Accession No.	Swiss-Prot: Q10589NCBI Gene ID: 684
Uniprot	Q10589
GeneID	684;
SDS-PAGE MW	20kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of SKOV3 cells using BST-2 Polyclonal Antibody

Note: This product is for in vitro research use only