

CD244 Polyclonal Antibody

Catalog No: #41980

Package Size: #41980-1 50ul #41980-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD244 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	CD244 Polyclonal Antibody detects endogenous levels of CD244 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CD244.
Target Name	CD244
Other Names	CD244; 2B4; Natural killer cell receptor 2B4; NK cell activation-inducing ligand; NAIL; NK cell type I receptor protein 2B4; NKR2B4; h2B4; CD244
Accession No.	Swiss-Prot: Q9BZW8NCBI Gene ID: 51744
Uniprot	Q9BZW8
GeneID	51744;
SDS-PAGE MW	42kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.


Not yet tested in other applications.

Images



Western Blot analysis of KB cells using CD244 Polyclonal Antibody

Western Blot analysis of KB cells using CD244 Polyclonal Antibody



Note: This product is for in vitro research use only