

IL-13Rα2 Polyclonal Antibody

Catalog No: #41981



Package Size: #41981-1 50ul #41981-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	IL-13Rα2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	IL-13Rα2 Polyclonal Antibody detects endogenous levels of IL-13Rα2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human IL-13Rα2.
Target Name	IL-13Rα2
Other Names	IL13RA2; IL13R; Interleukin-13 receptor subunit alpha-2; IL-13 receptor subunit alpha-2; IL-13R subunit alpha-2; IL-13R-alpha-2; IL-13RA2; Interleukin-13-binding protein; CD213a2
Accession No.	Swiss-Prot: Q14627NCBI Gene ID: 3598
Uniprot	Q14627
GeneID	3598;
SDS-PAGE MW	44kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.


Not yet tested in other applications.

Images



Western Blot analysis of HEPG2-UV cells using IL-13Rα2 Polyclonal Antibody

Western Blot analysis of HepG2(UV) L929 cells using IL-13R α 2 Polyclonal Antibody



Note: This product is for in vitro research use only