

CD305 Polyclonal Antibody

Catalog No: #41998

Package Size: #41998-1 50ul #41998-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD305 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	CD305 Polyclonal Antibody detects endogenous levels of CD305 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CD305.
Target Name	CD305
Other Names	LAIR1; CD305; Leukocyte-associated immunoglobulin-like receptor 1; LAIR-1; hLAIR1; CD305
Accession No.	Swiss-Prot: Q6GTX8NCBI Gene ID: 3903
Uniprot	Q6GTX8
GeneID	3903;
SDS-PAGE MW	32kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of KB cells using CD305 Polyclonal Antibody



Western Blot analysis of KB NIH-3T3 cells using CD305 Polyclonal Antibody

Note: This product is for in vitro research use only