

Neuregulin-3 Polyclonal Antibody

Catalog No: #42012

Package Size: #42012-1 50ul #42012-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neuregulin-3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms
Specificity	Neuregulin-3 Polyclonal Antibody detects endogenous levels of Neuregulin-3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Neuregulin-3.
Target Name	Neuregulin-3
Other Names	NRG3; Pro-neuregulin-3, membrane-bound isoform; Pro-NRG3
Accession No.	Swiss-Prot: P56975NCBI Gene ID: 10718
Uniprot	P56975
GeneID	10718;
SDS-PAGE MW	80kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.


Not yet tested in other applications.

Images



Western Blot analysis of MOUSE-BRAIN cells using
Neuregulin-3 Polyclonal Antibody

Western Blot analysis of mouse brain K562 cells using
Neuregulin-3 Polyclonal Antibody



Note: This product is for in vitro research use only