

# Complement C1q tumor necrosis factor-related protein 3 Polyclonal Antibody

Catalog No: #42089

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Complement C1q tumor necrosis factor-related protein 3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total Complement C1q tumor necrosis factor-related protein 3 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human Complement C1q tumor necrosis factor-related protein 3 protein
Target Name	Complement C1q tumor necrosis factor-related protein 3
Other Names	Collagenous repeat-containing sequence 26 kDa protein, Secretory protein CORS26, C1QTNF3, CTRP3
Accession No.	Swiss-Prot#: Q9BXJ4
Uniprot	Q9BXJ4
GeneID	114899;
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

## Application Details

Immunohistochemistry: 1:20 - 1:200

## Images



Immunohistochemical analysis of paraffin-embedded human colorectal carcinoma using #42089 at dilution of 1:20.

## Background

C1QTNF3 is a secretory protein produced by chondrogenic precursors & proliferating chondrocytes, and belongs to a novel C1q family of proteins.

## References

---

[1] "Homo sapiens complement-c1q tumor necrosis factor-related protein."Piddington C.S., Bishop P.Submitted (DEC-2000) [2] "Molecular cloning, chromosomal localization, and genomic structure of the human CORS26 gene."Maeda T., Hayashi A., Saito T.Submit

---

Note: This product is for in vitro research use only