

Clathrin heavy chain 2 Polyclonal Antibody

Catalog No: #42127

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Clathrin heavy chain 2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total Clathrin heavy chain 2 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human Clathrin heavy chain 2 protein
Target Name	Clathrin heavy chain 2
Other Names	Clathrin heavy chain on chromosome 22, CLH-22, CLTCL1, CLH22, CLTCL, CLTD
Accession No.	Swiss-Prot#: P53675
Uniprot	P53675
GeneID	8218;
Calculated MW	187kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Western blotting: □1:500 - 1:1000

Images



All lanes: Clathrin heavy chain 2 antibody at 2ug/ml

Lane 1: Jurkat whole cell lysate

Lane 2: NIH3T3 whole cell lysate

Lane 3: A431 whole cell lysate

secondary

Goat polyclonal to rabbit at 1/10000 dilution

predicted band size : 187kDa

observed band size : 187kDa

Background

Clathrin is the major protein of the polyhedral coat of coated pits and vesicles. Two different adapter protein complexes link the clathrin lattice either to the plasma membrane or to the trans-Golgi network

References

[1]"Quantitative phosphoproteomics reveals widespread full phosphorylation site occupancy during mitosis."Olsen J.V., Vermeulen M., Santamaria A., Kumar C., Miller M.L., Jensen L.J., Gnad F., Cox J., Jensen T.S., Nigg E.A., Brunak S., Mann M.S

Note: This product is for in vitro research use only