

60S ribosomal protein L8 Polyclonal Antibody

Catalog No: #42136

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	60S ribosomal protein L8 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total 60S ribosomal protein L8 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human 60S ribosomal protein L8 protein
Target Name	60S ribosomal protein L8
Other Names	RPL8
Accession No.	Swiss-Prot#: P62917
Uniprot	P62917
GeneID	6132;
Calculated MW	28kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500 - 1:1000

Images



All lanes : 60S ribosomal protein L8 antibody at 2ug/ml
 Lane 1 : EC109 whole cell lysate
 Lane 2 : 293T whole cell lysate
 Predicted band size : 28 kDa
 Observed band size : 28 kDa

Background

This protein can be partially incorporated into E.coli polysomes in vivo, indicating it can replace the endogenous protein.

References

[1]"The status, quality, and expansion of the NIH full-length cDNA project: the Mammalian Gene Collection (MGC)."The MGC Project TeamGenome Res. 14:2121-2127(2004) [2]"A map of 75 human ribosomal protein genes." Kenmochi N., Kawaguchi T., Rozen S., Davi

Note: This product is for in vitro research use only