

DTW domain-containing protein 1 Polyclonal Antibody

Catalog No: #42153

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DTW domain-containing protein 1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total DTW domain-containing protein 1 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant Human DTW domain-containing protein 1 protein
Target Name	DTW domain-containing protein 1
Other Names	DTWD1, MDS009
Accession No.	Swiss-Prot#: Q8N5C7
Uniprot	Q8N5C7
GeneID	56986;
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:20 - 1:200

Images



Immunohistochemical analysis of paraffin-embedded mouse liver using #42153 at dilution of 1:100 .

References

- [1] Novel genes expressed in hematopoietic stem/progenitor cells from myelodysplastic syndromes patient. Huang C., Zhang C., Tu Y., Gu W., Wang Y., Han Z., Chen Z., Zhou J., Gu J., Huang Q., Yu Y., Xu S., Ren S., Fu G. Submitted (JUL-1999) [2] Analysis of

Note: This product is for in vitro research use only