

Guanidinoacetate N-methyltransferase Polyclonal Antibody

Catalog No: #42175

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Guanidinoacetate N-methyltransferase Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB
Species Reactivity	Rt
Specificity	The antibody detects endogenous level of total Guanidinoacetate N-methyltransferase polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant rat Guanidinoacetate N-methyltransferase protein
Target Name	Guanidinoacetate N-methyltransferase
Other Names	Gamt
Accession No.	Swiss-Prot#: P10868
Uniprot	P10868
Calculated MW	26kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500 - 1:1000

Images



All lanes: Guanidinoacetate N-methyltransferase antibody at 2ug/ml+293T
 Secondary
 Goat polyclonal to Rabbit IgG at 1/10000 dilution
 Predicted band size: 26kDa
 Observed band size: 26kDa

References

[1] Molecular cloning, sequence analysis, and expression in Escherichia coli of the cDNA for guanidinoacetate methyltransferase from rat liver. Ogawa H., Date T., Gomi T., Konishi K., Pitot H.C., Cantoni G.L., Fujioka M. Proc. Natl. Acad. Sci. U.S.A. 85:694-

Note: This product is for in vitro research use only