

Probable histidine--tRNA ligase, mitochondrial Polyclonal Antibody

Catalog No: #42199

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Probable histidine--tRNA ligase, mitochondrial Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous level of total Probable histidine--tRNA ligase, mitochondrial polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human Probable histidine--tRNA ligase, mitochondrial protein
Target Name	Probable histidine--tRNA ligase, mitochondrial
Other Names	Histidine--tRNA ligase-like Histidyl-tRNA synthetase HARS2 HARSL, HARSR, HO3
Accession No.	Swiss-Prot#: P49590
Uniprot	P49590
GeneID	23438;
Calculated MW	57kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Western blotting: □1:500 - 1:1000

Images



All lanes : Probable histidine--tRNA ligase, mitochondrial antibody at 2ug/ml+mouse kidney tissue
 Secondary
 Goat polyclonal to Rabbit IgG at 1/10000 dilution
 Predicted band size:57kDa
 Observed band size:57kDa

References

[1]A novel gene oriented in a head-to-head configuration with the human histidyl-tRNA synthetase (HRS) gene encodes an mRNA that predicts a

Note: This product is for in vitro research use only