

Urokinase plasminogen activator surface receptor Polyclonal Antibody

Catalog No: #42294

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Urokinase plasminogen activator surface receptor Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB
Species Reactivity	Hu Rt
Specificity	The antibody detects endogenous level of total Urokinase plasminogen activator surface receptor polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant rat Urokinase plasminogen activator surface receptor ptoein
Target Name	Urokinase plasminogen activator surface receptor
Other Names	CD_antigen=CD87 Plaur
Accession No.	Swiss-Prot#: P29598
Uniprot	P29598
Calculated MW	34kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500 - 1:1000

Images



All lanes:Urokinase plasminogen activator surface receptor antibody at 2ug/ml
 Lane 1:293T whole cell lysate
 Lane 2:rat brain tissue
 secondary
 Goat polyclonal to rabbit at 1/10000 dilution
 predicted band size :34kDa
 observed band size :34kDa

Background

Specifically cleaves the zymogen plasminogen to form the active enzyme plasmin.

References

[1] Isolation and characterization of multiple isoforms of the rat urokinase receptor in osteoblasts. Rabbani S.A., Rajwans N., Achbarou A., Murthy K.K., Goltzman D. FEBS Lett. 338:69-74 (1994) [2] DeYoung M.B., Wohlgemuth J., Dichek D.A. Submitted (JUN-1997)

Note: This product is for in vitro research use only