

40S ribosomal protein S11 Polyclonal Antibody

Catalog No: #42299

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	40S ribosomal protein S11 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total 40S ribosomal protein S11 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human 40S ribosomal protein S11 protein
Target Name	40S ribosomal protein S11
Other Names	RPS11
Accession No.	Swiss-Prot#: P62280
Uniprot	P62280
GeneID	6205;
Calculated MW	18kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details


Western blotting: □ 1:500 - 1:1000

Immunohistochemistry: 1:20 - 1:200

Images



All lanes: 40S ribosomal protein S11 antibody at 2ug/ml+293T whole cell lysate
Secondary
Goat polyclonal to Rabbit IgG at 1/10000 dilution
Predicted band size:18kDa
Observed band size:18kDa



Immunohistochemical analysis of paraffin-embedded human placenta using #42299 at dilution of 1:100.

References

[1] Sequence of a cloned cDNA encoding human ribosomal protein S11. Lott J.B., Mackie G.A. Nucleic Acids Res. 16:1205-1205(1988) [2] Gene organization and sequence of the region containing the ribosomal protein genes RPL13A and RPS11 in the human genome

Note: This product is for in vitro research use only