

Non-secretory ribonuclease Polyclonal Antibody

Catalog No: #42307

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Non-secretory ribonuclease Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total Non-secretory ribonuclease polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human Non-secretory ribonuclease protein
Target Name	Non-secretory ribonuclease
Other Names	Eosinophil-derived neurotoxin, RNase Upl-2, Ribonuclease 2, Ribonuclease US, RNASE2, EDN, RNS2
Accession No.	Swiss-Prot#: P10153
Uniprot	P10153
GeneID	6036;
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:20 - 1:200

Images



Immunohistochemical analysis of paraffin-embedded human kidney using #42307 at dilution of 1:100.

Background

This is a non-secretory ribonuclease. It is a pyrimidine specific nuclease with a slight preference for U. Cytotoxin and helminthotoxin. Selectively chemotactic for dendritic cells. Possesses a wide variety of biological activities.

References

[1] "Sequence of human eosinophil-derived neurotoxin cDNA: identity of deduced amino acid sequence with human nonsecretory ribonucleases." Hamann K.J., Barker R.L., Loegering D.A., Pease L.R., Gleich G.J. *Gene* 83:161-167(1989) [2] "Eosinophil cationic pro

Note: This product is for in vitro research use only