

C-C motif chemokine 24 Polyclonal Antibody

Catalog No: #42444

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	C-C motif chemokine 24 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total C-C motif chemokine 24 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human C-C motif chemokine 24 protein
Target Name	C-C motif chemokine 24
Other Names	MPIF2, SCYA24, CCL24, CK-beta-6, Eosinophil chemotactic protein 2, Eotaxin-2, Myeloid progenitor inhibitory factor 2, MPIF-2, Small-inducible cytokine A24
Accession No.	Swiss-Prot#: O00175
Uniprot	O00175
GeneID	6369;
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:20 - 1:200

Images



Immunohistochemical analysis of paraffin-embedded human lymph gland using #42444 at dilution of 1:100.

Background

Chemotactic for resting T-lymphocytes, and eosinophils. Has lower chemotactic activity for neutrophils but none for monocytes and activated lymphocytes. Is a strong suppressor of colony formation by a multipotential hematopoietic progenitor cell line. Binds to CCR3.

References

[1] "Molecular and functional characterization of two novel human C-C chemokines as inhibitors of two distinct classes of myeloid progenitors." Patel V.P., Kreider B.L., Li Y., Li H., Leung K., Salcedo T., Nardelli B., Pippalla V., Gentz S., Thotakura R.,

Note: This product is for in vitro research use only