

Golgi membrane protein 1 Polyclonal Antibody

Catalog No: #42527

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Golgi membrane protein 1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total Golgi membrane protein 1 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant Golgi membrane protein 1
Target Name	Golgi membrane protein 1
Other Names	PSEC0242, UNQ686/PRO1326, C9orf155, GOLPH2, GOLM1, Golgi membrane protein GP73, Golgi phosphoprotein 2
Accession No.	Swiss-Prot#: Q8NBJ4
Uniprot	Q8NBJ4
GeneID	51280;
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:20 - 1:200

Images



Immunohistochemical analysis of paraffin-embedded human prostate using #42527 at dilution of 1:100.

Background

Cellular response protein to viral infection.

References

[1] "GP73, a novel Golgi-localized protein upregulated by viral infection."Kladney R.D., Bulla G.A., Guo L., Mason A.L., Tollefson A.E., Simon D.J., Koutoubi Z., Fimmel C.J.Gene 249:53-65(2000) [PubMed: 10831838] [Abstract]Cited for: NUCLEOTIDE SEQUENCE [

Note: This product is for in vitro research use only