

STYK1 antibody

Catalog No: #22780

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	STYK1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 214 and 400 of Human STYK1
Target Name	STYK1
Accession No.	NCBI Gene ID: 55359NCBI mRNA#: NM_018423NCBI Protein#: NP_060893
Uniprot	Q6J9G0
GeneID	55359;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 48kd

Western blotting: 1:500-1:3000

Immunohistochemistry: 1:100-1:250

Images




Sample(30 ug whole cell lysate)

A: A431

B: MOLT4

7.5% SDS PAGE

Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded Breast CA, using STYK1 antibody at 1: 100 dilution.

Background

Receptor protein tyrosine kinases, like STYK1, play important roles in diverse cellular and developmental processes, such as cell proliferation, differentiation, and survival (Liu et al., 2004 [PubMed 15150103]).[supplied by OMIM]

Note: This product is for in vitro research use only