

ITPK1 antibody

Catalog No: #22783



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


Description

Product Name	ITPK1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 171 and 405 of ITPK1
Target Name	ITPK1
Accession No.	Swiss-Prot:Q13572Gene ID:3705
Uniprot	Q13572
GeneID	3705;
Concentration	0.5mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 46kd
Western blotting: 1:500-1:3000
Immunofluorescence: 1:100-1:200


Images



Sample (30 ug of whole cell lysate)  
A: A431  
10% SDS PAGE  
Primary antibody diluted at 1: 1000



Merged with DNA probe



Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using ITPK1 antibody at 1: 200 dilution.

---

Note: This product is for in vitro research use only