

Probable sarcosine oxidase Polyclonal Antibody

Catalog No: #42621

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Probable sarcosine oxidase Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB
Species Reactivity	Arabidopsis thaliana
Specificity	The antibody detects endogenous level of total Probable sarcosine oxidase polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant Arabidopsis thaliana Probable sarcosine oxidase protein
Target Name	Probable sarcosine oxidase
Other Names	At2g24580, F25P17.12
Accession No.	Swiss-Prot#: Q9SJA7
Uniprot	Q9SJA7
GeneID	816994;
Calculated MW	46kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Western blotting: □1:500 - 1:1000

Images



All lanes : Probable sarcosine oxidase antibody at 2ug/ml

Lane 1 : 293T whole cell lysate

Lane 2 : EC109 whole cell lysate

Secondary

Goat polyclonal to Rabbit IgG at 1/10000 dilution

Predicted band size :46KDa

Observed band size:46KDa

References

[1] Sequence and analysis of chromosome 2 of the plant Arabidopsis thaliana. Lin X., Kaul S., Rounsley S.D., Shea T.P., Benito M.-I., Town C.D., Fujii C.Y., Mason T.M., Bowman C.L., Barnstead M.E., Feldblyum T.V., Buell C.R., Ketchum K.A., Lee J.J., Ronnin

Note: This product is for in vitro research use only