

Rat Tumor Necrosis Factor Receptor Superfamily, Member 5 (TNFRSF5) ELISA Kit

Catalog No: #EK17018

Package Size: #EK17018-1 48T #EK17018-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Tumor Necrosis Factor Receptor Superfamily, Member 5 (TNFRSF5) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rattus norvegicus Rat
Other Names	CD40; CDW40; Bp50; P50; B-cell surface antigen CD40; CD40L receptor
Accession No.	UniProtKB:Q4QQW2Gene ID:171369
Uniprot	Q4QQW2
GeneID	171369;
Storage	12 months

Application Details

Detect Range:15.6-1000pg/mL

Sensitivity:5.2pg/mL

Sample Type:serum, plasma or other biological fluids.

Sample Volume: 100 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only