

## Mouse IgG

Catalog No: #IG102



Package Size: #IG102&lt;br /&gt;-1 100mg #IG102&lt;br /&gt;-2 200mg #IG102&lt;br /&gt;-3 1g

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse IgG
Purification	Affinity chromatography or caprylic acid-ammonium sulfate,85%
Conjugates	None
Concentration	20 mg/ml
Storage	-20°C

## Product Description

Source: Normal, healthy, unimmunized mouse serumAppearance: Liquid Filtration: 0.2um Buffer :10mM Phosphate buffered salinePBSpH7.40.Preservative: NoneStability: Stable for 24 monthsPrecautions: This product can be kept at 2-8°C for maximum 30 days. Avoid repeated freezing and thawing. If multiple usages are needed, it is recommended to be stored separately in smaller volumes.

Note: This product is for in vitro research use only