

PDE4C antibody

Catalog No: #22806

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PDE4C antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 325 and 637 of Human PDE4C
Target Name	PDE4C
Accession No.	Swiss-Prot:Q08493Gene ID:5143
Uniprot	Q08493
GeneID	5143;
Concentration	0.7mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 80kd


Western blotting: 1:500-1:3000

Immunohistochemistry: 1:100-1:250

Images



Sample(30 ug of whole cell lysate)
A: 293T
7.5% SDS PAGE
Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded endometrial, using PDE4C antibody at 1: 100 dilution.

Note: This product is for in vitro research use only