

Annexin A1 Antibody

Catalog No: #P1024

Package Size: #P1024-1 0.1ml #P1024-2 1ml

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Annexin A1 Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	IHC512
Applications	IHC
Species Reactivity	Human
Formulation	Tris Buffer, pH 7.3 - 7.7, with 1% BSA and
Storage	Store at 2-8C. Do not freeze.

Application Details

Recommended

working dilution range:1:100 - 1:200

Images



Annexin A1 Antibody (P1024) on Spleen

Product Description

Annexin A1 (ANXA1) is a membrane protein that plays a role in innate and adaptive immunity by controlling the biosynthesis of inflammation, prostaglandins, and leukotriene mediators. This target is overexpressed in 97% of all samples from patients with hairy cell leukemia, and is absent in other B-cell lymphomas. High ANXA1 expression is frequently associated with advanced stage esophageal and esophagogastric junction adenocarcinoma, and is also linked to advanced and metastatic disease states.

References

Tiaci E, et al. N Eng J Med. 2011; 364:2305-2315.

Wang K, et al. Clin Cancer Res. 2006; 12:4598-4604.

Note: This product is for in vitro research use only