

## Recombinant Human ANXA4

Catalog No: #GP10165



Package Size: #GP10165-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ANXA4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 179-314 amino acids of human annexin A4
Target Name	annexin A4
Other Names	ANX4; PIG28; ZAP36
Accession No.	Swissprot:P09525Gene Accession:BC000182
Uniprot	P09525
GeneID	307;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Annexin IV (ANX4) belongs to the annexin family of calcium-dependent phospholipid binding proteins. Although their functions are still not clearly defined, several members of the annexin family have been implicated in membrane-related events along exocytotic and endocytotic pathways. ANX4 has 45 to 59% identity with other members of its family and shares a similar size and exon-intron organization. Isolated from human placenta, ANX4 encodes a protein that has possible interactions with ATP, and has in vitro anticoagulant activity and also inhibits phospholipase A2 activity. ANX4 is almost exclusively expressed in epithelial cells.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only