

MYBPHL Antibody

Catalog No: #42913

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MYBPHL Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total MYBPHL protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human MYBPHL
Target Name	MYBPHL
Other Names	MBPHL; MYBPHL; myosin binding protein H-like; Myosin-binding protein H-like Gene Symbols: MYBPHL
Accession No.	Swiss-Prot#: A2RUH7Gene ID: 343263
Uniprot	A2RUH7
GeneID	343263;
Concentration	0.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human breast cancer tissue using #42913 at dilution 1/20.



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #42913 at dilution 1/20.

Background

This gene encodes a protein with two immunoglobulin superfamily domains and a fibronectin 3 domain. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only