

PRSS38 Antibody

Catalog No: #42924

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PRSS38 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total PRSS38 protein.
Immunogen Type	protein
Immunogen Description	Full length fusion protein of human PRSS38
Target Name	PRSS38
Other Names	MPN2
Accession No.	Swiss-Prot#: A1L453Gene ID: 339501
Uniprot	A1L453
GeneID	339501;
Concentration	1.3mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #42924 at dilution 1/30.



Immunohistochemical analysis of paraffin-embedded Human cervical cancer tissue using #42924 at dilution 1/30.

Background

The function of this protein remains unknown.

Note: This product is for in vitro research use only