

# Recombinant Human PGRMC2

Catalog No: #GP11861

Package Size: #GP11861-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human PGRMC2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 123 amino acids of human PGRMC2
Target Name	progesterone receptor membrane component 2
Other Names	DG6; PMBP
Accession No.	Swissprot:O15173Gene Accession:BC016692
Uniprot	O15173
GeneID	10424;
Storage	-20~-80°C, pH 7.6 PBS

## Background

PGRMC2 (progesterone receptor membrane component 2), also known as DG6 (steroid receptor protein DG6) or PMBP (progesterone membrane-binding protein), is a single pass membrane protein belonging to the cytochrome b5 family (MAPR (membrane associated progesterone receptor) subfamily). Expressed in sperm, PGRMC2 is believed to function as a steroid receptor and may participate in the progesterone-dependent sperm acrosome reaction.?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only