

## SND1 Antibody

Catalog No: #43043

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	SND1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total SND1 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human SND1
Target Name	SND1
Other Names	p100; TDRD11
Accession No.	Swiss-Prot#: Q7KZF4 Gene ID: 27044
Uniprot	Q7KZF4
GeneID	27044;
Calculated MW	102kd
Concentration	1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-1:5000

Immunohistochemistry: 1:25-1:100

## Images



Gel: 6%SDS-PAGE


Lysate: 40 µg

Lane 1-6: HeLa and Jurkat cell, mouse liver tissue, hepg2, A172 and NIH/3T3 cell

Primary antibody: 1/400 dilution

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 seconds



Immunohistochemical analysis of paraffin-embedded Human breast cancer tissue using #43043 at dilution 1/20. (Original magnification: B1F• 200)

## Background

Functions as a bridging factor between STAT6 and the basal transcription factor. Plays a role in PIM1 regulation of MYB activity. Functions as a transcriptional coactivator for the Epstein-Barr virus nuclear antigen 2 (EBNA2).

Note: This product is for in vitro research use only