

SKI Antibody

Catalog No: #43088

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SKI Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total SKI protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human SKI
Target Name	SKI
Other Names	SGS; SKV
Accession No.	Swiss-Prot#: P12755Gene ID: 6497
Uniprot	P12755
GeneID	6497;
Concentration	2.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #43088 at dilution 1/40.



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #43088 at dilution 1/40.

Background

This gene encodes the nuclear protooncogene protein homolog of avian sarcoma viral (v-ski) oncogene. It functions as a repressor of TGF-beta signaling, and may play a role in neural tube development and muscle differentiation.?

Note: This product is for in vitro research use only