

Recombinant Human WDFY1

Catalog No: #GP12541

Package Size: #GP12541-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human WDFY1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 211-410 amino acids of human WDFY1
Target Name	WD repeat and FYVE domain containing 1
Other Names	WDF1; FENS-1; ZFYVE17
Accession No.	Swissprot:Q8IWB7Gene Accession:BC040525
Uniprot	Q8IWB7
GeneID	57590;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a phosphatidylinositol 3-phosphate binding protein, which contains a FYVE zinc finger domain and multiple WD-40 repeat domains. When exogenously expressed, it localizes to early endosomes. Mutagenesis analysis demonstrates that this endosomal localization is mediated by the FYVE domain.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only