

SAH3 antibody

Catalog No: #22835



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


Description

Product Name	SAH3 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 297 and 547 of Human AHCYL2
Target Name	SAH3
Accession No.	Swiss-Prot:Q96HN2Gene ID:23382
Uniprot	Q96HN2
GeneID	23382;
Concentration	0.7mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 67kd  
Western blotting: 1:500-1:3000  
Immunohistochemistry: 1:100-1:250

Images



Sample (30 ug of whole cell lysate)  
A: Raji  
7.5% SDS PAGE  
Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded BT474 xenograft, using AHCYL2 antibody at 1: 100 dilution.

---

Note: This product is for in vitro research use only