

Recombinant Human CHAC2

Catalog No: #GP12765

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CHAC2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	ChaC cation transport regulator homolog 2
Accession No.	Swissprot:Q8WUX2Gene Accession:BC019239
Uniprot	Q8WUX2
GeneID	494143;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a gamma-glutamyl cyclotransferase that catalyzes the conversion of glutathione to 5-oxoproline and cysteinylglycine. It is thought that this gene is upregulated in response to endoplasmic reticulum stress and that the glutathione depletion enhances apoptosis.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only