

# Recombinant Human PSD4

Catalog No: #GP13212

Package Size: #GP13212-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human PSD4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 857-1056 amino acids of human PSD4
Target Name	pleckstrin and Sec7 domain containing 4
Other Names	TIC; EFA6B
Accession No.	Swissprot:Q8NDX1Gene Accession:BC035307
Uniprot	Q8NDX1
GeneID	23550;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Guanine nucleotide exchange factor for ARF6 and ARL14/ARF7. Through ARL14 activation, controls the movement of MHC class II-containing vesicles along the actin cytoskeleton in dendritic cells. Involved in membrane recycling. Interacts with several phosphatidylinositol phosphate species, including phosphatidylinositol 3,4-bisphosphate, phosphatidylinositol 3,5-bisphosphate and phosphatidylinositol 4,5-bisphosphate.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only