

AGPAT6 Antibody

Catalog No: #43170

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	AGPAT6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total AGPAT6 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human AGPAT6
Target Name	AGPAT6
Other Names	LPAATZ; TSARG7; 1-AGPAT 6; LPAAT-zeta
Accession No.	Swiss-Prot#: Q86UL3Gene ID: 137964
Uniprot	Q86UL3
GeneID	137964;
Calculated MW	52kd
Concentration	1.6mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:200-1:1000

Immunohistochemistry: 1:25-1:100

Images



Gel: 8%SDS-PAGE
Lysate: 40 µg
Lane: Mouse plasma tissue
Primary antibody: 1/400 dilution
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution
Exposure time: 5 minutes

Background

Lysophosphatidic acid acyltransferases (EC 2.3.1.51) catalyze the conversion of lysophosphatidic acid (LPA) to phosphatidic acid (PA). LPA and PA

are involved in signal transduction and lipid biosynthesis.

Note: This product is for in vitro research use only