

## TSPYL1 Antibody

Catalog No: #43171

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	TSPYL1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total TSPYL1 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human TSPYL1
Target Name	TSPYL1
Other Names	TSPYL
Accession No.	Swiss-Prot#: Q9H0U9Gene ID: 7259
Uniprot	Q9H0U9
GeneID	7259;
Calculated MW	49kd
Concentration	1.2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:200-1:1000

Immunohistochemistry: 1:25-1:100

## Images



Gel: 8%SDS-PAGE  
Lysate: 40 µg  
Lane: Mouse brain tissue  
Primary antibody: 1/200 dilution  
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution  
Exposure time: 2 minutes

## Background

The protein encoded by this gene is found in the nucleolus and is similar to that of a family of genes on the Y-chromosome. This gene is intronless.

Defects in this gene are a cause of sudden infant death with dysgenesis of the testes syndrome (SIDDT).?

---

Note: This product is for in vitro research use only